



**RETURN PRESCRIPTION TO**  
9129 Lurline Ave.  
Chatsworth, CA 91311  
800-423-3270  
www.SMLglobal.com



**PLEASE SEND:**

- ☐ MAILING MATERIALS  
☐ PRODUCTS & SUPPLIES  
☐ DIAGNOSTIC SERVICES  
☐ C.E. COURSES



**ADDITIONAL SERVICES\***

- ☐ RETURN DUPLICATE SET OF MODELS  
☐ APPLIANCE INSURANCE

ACCOUNT #

DOCTOR NAME

OFFICE ADDRESS

CITY

STATE

ZIP CODE

PATIENT'S FIRST NAME

PATIENT'S LAST NAME

AGE

OFFICE EMAIL ADDRESS

DATE OF BIRTH

OFFICE PHONE NUMBER

**DUE DATE** – MUST BE A MINIMUM OF ONE DAY PRIOR TO YOUR PATIENT'S APPOINTMENT

(LAB USE ONLY)

- ☐ EMERGENCY SERVICE FOR APPLIANCES\*  
(24 to 48 Hrs. Processing)  
☐ PATIENT WILL BE APPOINTED AFTER APPLIANCE ARRIVES

☐ S.I.

**INDIRECT BONDING SYSTEM - THE COMPLETE PACKAGE**  
Includes: Digital Study Models with IPR Analysis and Individual Patient Kit (see reverse for details), plus Indirect Bonding Tray(s)  
☐ Composite Brackets ☐ Stainless Steel Brackets  
☐ Upper AND Lower ☐ Upper ☐ Lower  
☐ UD 018 Bracket ☐ UR2 ☐ UL2  
☐ Bracket for Shallow/Open Bite



**INDIRECT BONDING SYSTEM - INDIVIDUAL REQUEST OPTIONS**  
Indirect Bonding Tray(s) ☐ Composite Brackets ☐ Stainless Steel Brackets  
☐ Upper AND Lower ☐ Upper ☐ Lower  
☐ Digital Study Models with IPR Analysis  
☐ Individual Patient Kit (see reverse for details)  
☐ QCB Ceph Analysis

**APPLIANCES\***

- ☐ QCB Splint (includes 3x3 and Zendura Retainer) ☐ Hawley  
☐ 3x3 ☐ Myofunctional Tongue Bead  
☐ Zendura Retainer ☐ Nitanium Molar Rotator  
☐ Inman Aligners ☐ Nitanium Palatal Expander  
☐ Spring Retainer  
\* FEES APPLY

**SPECIAL INSTRUCTIONS:**

**UPPER** # \_\_\_\_\_ **LOWER** # \_\_\_\_\_  
☐ FIXED from Principles of Appliance Therapy textbook  
☐ REMOVABLE  
☐ FIXED from Principles of Appliance Therapy textbook  
☐ REMOVABLE



**BRACKET HEIGHT PRESCRIPTION**

|                                  |       |       |       |       |       |       |       |       |       |       |       |       |                                  |
|----------------------------------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|----------------------------------|
| Custom Height<br>Standard Height | 3.5mm | 4.0mm | 4.5mm | 5.0mm | 4.0mm | 4.5mm | 4.5mm | 4.0mm | 5.0mm | 4.5mm | 4.0mm | 3.5mm | Custom Height<br>Standard Height |
| R                                |       |       |       |       |       |       |       |       |       |       |       |       | L                                |
| Standard Height<br>Custom Height | 3.5mm | 4.0mm | 4.5mm | 5.0mm | 4.0mm | 4.0mm | 4.0mm | 4.0mm | 5.0mm | 4.5mm | 4.0mm | 3.5mm | Standard Height<br>Custom Height |

☐ Additional Instructions Attached

SIGNATURE

LICENSE NUMBER

**GO GREEN! PLEASE SCAN OR MAKE A COPY OF THIS PRESCRIPTION FORM FOR YOUR RECORDS**

**BEFORE SUBMITTING TO LAB:**

- ☐ **PRESCRIPTION** - Make sure all appropriate sections are completed.  
☐ **STONE MODELS** - Be sure to get doctor's final approval on models (to ensure accuracy and completeness). Trim models as small as possible.  
☐ **DIGITAL RECORDS** - If applicable, send digital patient files to www.SMLglobal.com/digital

- ☐ **ACCURATE CONSTRUCTION BITE** - Include for all cases where acrylic occlusal coverage or mandibular repositioning is required.  
☐ **PACKAGING** - Sturdy cardboard box (provided upon request) is required. Fill box completely with packing material. Wrap models carefully and individually.

**Terms and conditions on reverse**

# Quick Cosmetic Braces

## ☐ Individual Patient Kit

**Brackets:** Composite Brackets (1 set) - Upper and Lower

**Archwires:** One each of the following

|                                                 |                                                 |
|-------------------------------------------------|-------------------------------------------------|
| Upper .014 Tooth-Colored NiTi Wire              | Lower .014 Tooth-Colored NiTi Wire              |
| Upper .016 Tooth-Colored NiTi Wire              | Lower .016 Tooth-Colored NiTi Wire              |
| Upper .018 Tooth-Colored NiTi Wire              | Lower .018 Tooth-Colored NiTi Wire              |
| Upper .018 Tooth-Colored Stainless Steel Wires  | Lower .018 Tooth-Colored Stainless Steel Wires  |
| Upper 16x22 Tooth-Colored Stainless Steel Wires | Lower 16x22 Tooth-Colored Stainless Steel Wires |
| Upper 17x25 Tooth-Colored NiTi Wire             |                                                 |
| Upper 17x25 Tooth-Colored Stainless Steel Wires |                                                 |
| Upper 016 NiTi Reverse Curve - Straight Leg     |                                                 |

**Auxillaries:** Clear Ligature Ties (1 pack), Clear Chain Elastics (1 roll), Bracket Bonding Agent Ortho Solo (1 set-up), Bracket Adhesive Enlight (1 set-up), Diamond Disc w/ Clear Guard (1), Diamond Strip (1), Satin Diamond Bur (1), Shorty Twist Ties (pk 10), Patient Wax (1 each).

## Instruments and Supplies:

|                                                                   |                     |
|-------------------------------------------------------------------|---------------------|
| <input type="checkbox"/> 230-223 Mathieu Plier, Narrow Tip        | \$ 26.00 each       |
| <input type="checkbox"/> 230-620 Pin and Ligature Cutter          | \$ 72.65 each       |
| <input type="checkbox"/> 220-205 Distal End Cutter w/ safety hold | \$ 99.75 each       |
| <input type="checkbox"/> 220-319 Deluxe D.B. Bracket Holder       | \$215.00 Pack of 10 |
| <input type="checkbox"/> 230-760 Tube Placement Tool              | \$ 46.00 each       |
| <input type="checkbox"/> 220-310 Ligature Director                | \$ 18.40 each       |
| <input type="checkbox"/> 230-622 Bracket Height Gauge             | \$ 53.20 each       |
| <input type="checkbox"/> 220-312 Bracket Remover                  | \$ 80.50 each       |
| <input type="checkbox"/> 230-218 Weingart Plier                   | \$ 65.50 each       |

## Masters Series “Add-Ons”

|                                                                                                                     |                    |
|---------------------------------------------------------------------------------------------------------------------|--------------------|
| <input type="checkbox"/> 230-756 Step Bend Plier (1mm)*                                                             | \$ 79.95 each      |
| <input type="checkbox"/> 230-786 Hollow Chop Plier                                                                  | \$ 76.25 each      |
| <input type="checkbox"/> 420-128 Gruin Lock Wrench                                                                  | \$ 11.25 each      |
| <input type="checkbox"/> 420-127 Gurin Locks                                                                        | \$ 20.30 Pack of 2 |
| <input type="checkbox"/> 420-700 Clear Composite .022 Slot with Hooks and Buccal Tubes<br>(One upper and lower set) | \$ 85.00           |

---

---

### TERMS AND CONDITIONS

#### LABORATORY APPLIANCES

##### TERMS:

All invoices are due 15 days from invoice. At day 30, credit card on file will be charged. We accept Mastercard, Visa, American Express, and Discover. A 1.5% interest charge (18% per year) will be added to all invoices not paid by the due date. If legal action is required to obtain payment, SML is entitled to attorney fees.

##### LIABILITY RELEASE STATEMENT

SML provides appliances and laboratory services as prescribed by a licensed Dental Practitioner. We can assume no responsibly for techniques used and their use and/or misuse by the prescribing doctor, staff, or their patients.

##### APPLIANCE WARRANTY AND CONDITIONS

Our ability to provide a quality appliance begins with YOU. Please take the time to provide us an accurate set of working casts and a construction bite. Although we pride ourselves in our craftsmanship our appliances are only as good as the records provided for their fabrication.

SML can only guarantee that our custom made appliances will fit the working cast(s) provided for their construction. Materials and workmanship on all appliances are guaranteed for 90 days. If an appliance fails within this period, SML will remake or repair the appliance at no charge. This warranty does not cover appliance loss, patient abuse, or changes in the dentition during this period that would necessitate the need for a new appliance. All returns are subject to taxes, as well as regulatory compliance fees, model pour-up fees, and shipping fees.

##### SLEEP APPLIANCE WARRANTY AND CONDITIONS

Hardware and workmanship on all custom made sleep appliances are guaranteed for a period of 3 years. This warranty does not cover appliance loss, delamination of hard/soft material, patient abuse, or changes in the dentition during this period that would necessitate the need for a new appliance. All returns are subject to taxes, as well as regulatory compliance fees, model pour-up fees, and shipping fees.

##### REMAKES

**IMPORTANT NOTE:** While SML understands that many patients depend upon their devices for improved and continued health, requests for a total remake - while the patient continues to use the current appliance - should be neither expected by the dentist nor promised to the patient.

If an appliance does not fit your patient, follow these steps:

1. Send us a new accurate cast(s) or a digital impression along with a new bite registration.
2. Return the appliance that needs to be remade along with the original working casts used in its fabrication. These casts were returned to you with the original shipment of the appliance!
3. If the returned appliance does not fit the patient and does not fit the original working casts SML will fabricate a new appliance on your new casts at NO CHARGE.
4. If the appliance does not fit the patient but does fit the returned original working cast, SML will fabricate

a new appliance on your new casts and charges will be incurred at our usual and customary fees.

5. Occasionally the working casts used in fabrication may become damaged to a point that you will be unable to check the accuracy and fit of our work. Should this happen SML will document this occurrence and return a note with the appliance indicating “Damaged models during processing”. If the appliance does not fit the patient SML will remake the appliance at NO CHARGE. Simply follow steps 1 and 2 and write on the lab prescription slip that the case is a “broken cast remake”.

##### WHAT IS NOT COVERED BY WARRANTY

- Cash refund or a credit for custom dental device
- Cost incurred for patient canceling treatment
- Incidental or consequential damages, including inconvenience, lost wages, chair time, pain and suffering
- Improper adjustment of the appliance
- Concerns expressed to the doctor regarding impressions, bite registration, questionable indications and authorization for appliance manufacture
- Repairs and adjustments to reset a construction bite
- Repairs resulting from an accident, neglect, abuse or improper hygiene
- Changes in dentition; loss or removal of teeth; new restorations; failure of supportive tooth or tissue structures
- Partial or complete fabrication by a lab other than SML
- Expedited shipping costs
- De-lamination of hard/soft material
- Wear and tear from clenching/grinding

##### WHAT IS REQUIRED FOR WARRANTY COVERAGE

**NOTE:** Revamping or modifying of a broken device is considered standard procedure for warranty cases.

- Original appliance, impressions/models and bite must be returned to SML
- In the event of a fit problem, original impressions/models and bite must be returned in addition to new impressions and bite in the interest of ensuring accurate modifications and repairs to the original appliance.
- If new casts or impressions are required, simple refitting will be attempted under warranty. Should this necessitate a charge for unwarranted reasons, minimal or reasonable charges will be kept to a minimum or reasonably reduced level.
- Replacement screws and/or hinges may be provided in exchange for lost or defective screws and hinges - and will be credited upon return of the old parts.

##### WARNING

Many appliances are fabricated from stainless steel, nickel titanium and acrylic. Stainless steel contains small amounts of nickel and chromium. Nickel titanium contains nickel. Acrylic is processed with methyl methacrylate. A small number of the population is known to be allergic to these materials. Should an allergic reaction occur, advise the patient to consult a physician.